

Products for Biotechnology With Magnetic Porous Glass (MPG®)

Protocol No.: 47.4
Product: MPG® Goat anti-Mouse IgG (H&L) (10 mg/ml, 4-6 x 10⁷ particles/ml)
Procedure: Binding Mouse IgG
Storage: 2-8°C, Do Not Freeze.

PRODUCT NUMBER	DESCRIPTION	VOLUME
G-AB11020	MPG® Goat anti-Mouse IgG (H&L)	2 ml (20 mg)
G-AB11100		10 ml (100 mg)

MPG® Goat anti-Mouse IgG is a suspension of spheroidal totally porous borosilicate glass particles embedded with superparamagnetic iron oxide to which affinity purified Goat anti-Mouse IgG antibody has been covalently attached. This Goat secondary antibody binds both the heavy and light chains of all of the Mouse IgG subclasses (IgG₁, IgG_{2a}, IgG_{2b}, and IgG₃). Human cross reactivity is minimal. MPG® Secondary Antibody particles eliminate centrifugation, allow easier scale-up and automation, and give flexibility in assay configuration. These particles can be used as a secondary antibody in single or automated enzyme and radio immunoassays utilizing a Mouse IgG monoclonal or polyclonal antibody. They are highly recommended for separating unbound labeled tracer. These particles are excellent in small-scale immunoaffinity purifications or immunoprecipitation reactions for gel or other physical analyses. Use of magnetic particles not only reduces assay time, but as a result increases performance. These particles should be used in conjunction with a magnetic particle separator to achieve purification in minutes. The binding capacity of this antibody may vary among IgG isotypes and analyte bound IgG. In these instances, the researcher may have to determine actual binding capacities for their particular system. Check the accompanying certificate of analysis for lot specific information. The optimal dilution for specific applications should be determined by the researcher.

Store the vial at 2-8° C. This product contains **sodium azide** as a preservative. Sodium azide may react with lead and copper plumbing to form explosive metal azides. On disposal of material, flush with large volumes of water to prevent azide accumulation. Do not freeze or dry, as loss of binding activity and aggregation will occur. Dilute only prior to immediate use. Expiration date is six (6) months from date of opening product.

General Procedure

Materials: *MPG® Goat anti-Mouse IgG is supplied as a 1.0% suspension containing:*
Buffer: 0.01 M Phosphate, 0.15 M Sodium Chloride, pH 7.4
Stabilizer: 1 mg/ml Bovine Serum Albumin (BSA) – Immunoglobulin and protease free
Preservative: 0.02% (w/v) Sodium Azide

2 N Hydrochloric Acid (HCl)	Sodium Phosphate, Dibasic (Na ₂ HPO ₄)
Potassium Phosphate, Monobasic (KH ₂ PO ₄)	Tween 20
Potassium Chloride (KCl)	1.5 ml Microcentrifuge Tubes
Glycine (C ₂ H ₅ NO ₂)	Magnetic Particle Separator, Prod. No. MPS0301 or MPS0001
Sodium Chloride (NaCl)	Low Speed Rotator
Bovine Serum Albumin (BSA)	Pipettes and Pipette Tips
Deionized Water (dH ₂ O)	Vortex Mixer
Sodium Azide (NaN ₃)	Primary Antibody

Solution

1M Na₂HPO₄

Binding Buffer
(Phosphate-buffered saline pH 7.4,
{PBS})

Washing Buffer
(PBS, 0.05% Tween 20)

Storage Buffer*
(PBS, pH 7.4, 0.1% BSA, 0.02% NaN₃)

Elution Buffer
(0.1M Glycine-HCl, pH 2.5)

Preparation

Dissolve 142 g Na₂HPO₄ in 800 ml dH₂O and adjust volume to 1000 ml with dH₂O.

Dissolve 8 g NaCl, 0.2 g KCl, 1.17 g Na₂HPO₄ and 0.24 g KH₂PO₄ in 800 ml dH₂O. Adjust the pH to 7.4 with 2 N HCl and adjust volume to 1000 ml with dH₂O.

Mix 0.25 gm Tween 20 in 500 ml PBS.

Dissolve 100 mg BSA and 20 mg NaN₃ in 100 ml of Binding Buffer.

Dissolve 7.5 g Glycine in 500 ml dH₂O. Adjust pH to 2.5 with HCl and adjust volume to 1000 ml with H₂O.

**If the desired antibody cross-reacts with BSA, serum albumin of other species should be used in Storage Buffer.*

NOTE: These buffer formulations are only intended as a general guideline. Feel free to substitute and/or modify buffers as your application demands.

Preparation of MPG® Goat anti-Mouse IgG Particles:

Depending on your application a **preliminary wash step** using an appropriate buffer may be necessary. Washing will remove various additives including EDTA and anti-microbials, and concentrates the particles.

NOTE: For cell separations and further culturing, it is recommended that the particles be washed 5X in sterile medium prior to use.

1. Vortex the MPG® Goat anti-Mouse IgG to fully suspend the particles. Add the desired amount of particles to a 1.5 ml microcentrifuge tube. Magnetically separate the particles from the solution by placing the tube in a Magnetic Particle Separator for at least 30 seconds and carefully remove the supernatant by aspiration, with a pipette, while the tube remains in the particle separator.
2. Resuspend the MPG® Goat anti-Mouse IgG particles in Binding Buffer to original starting concentration (or more) and mix well. Magnetically separate and aspirate the supernatant. Repeat this step two more times.

Attachment of an Antibody to MPG® Goat anti-Mouse IgG Particles:

(The volumes and concentrations can be adjusted to suit particular applications.)

1. Resuspend the particles to give 0.5 mg solids/ml in a solution of Binding Buffer containing primary antibody. The antibody concentration should be determined empirically, but can be based on the binding capacity of the particles, as shown on the Certificate of Analysis.
2. Incubate at room temperature for 30 minutes with gentle mixing on low speed rotator.
3. Add Washing Buffer to the antibody-bound MPG® Goat anti-Mouse IgG particles to a concentration of 0.5 mg solids/ml and mix well. Magnetically separate and aspirate the supernatant. Repeat four more times. The antibody-bound MPG® Goat anti-Mouse IgG particles are now ready for use in your specific application.

Storage of Primary Antibody Coated Particles:

For storage, add 1 ml of Storage Buffer to the antibody-bound MPG® Goat anti-Mouse IgG particles and mix well. Magnetically separate and aspirate the supernatant. Resuspend the antibody-bound MPG® Goat anti-Mouse IgG particles in 1 ml of Storage Buffer to the desired storage concentration (often 0.5 mg/ml) and store at 4°C.

Elution of Primary Antibody from MPG® Goat anti-Mouse IgG Particles:

Break antibody/antigen interaction by suspending the complex in Elution Buffer for 5 minutes, with gentle mixing, at room temperature. Any convenient volume can be used, but the minimum should be 10-15 µl in order to adequately cover the particles and inside surface of the microfuge tube. Magnetic particles can then be separated and the supernatant, containing analyte, can be used in your application. Add 0.1 volume of 1 M Na₂HPO₄ to the recovered supernatant to neutralize the pH, if your analyte is acid labile.

**FOR TECHNICAL SERVICE ON THIS OR ANY OTHER PureBiotech PRODUCT CALL 866-252-7771
or e-mail us at info@purebiotechllc.com.**

For in vitro research use only.

MPG® is a registered trademark of Millipore Corporation. Magnetic Porous Glass and certain applications in which it is used are covered by U.S. Patent 5,601,979; 5,610,274 and 5,734,020 owned by Millipore Corporation and licensed by PureBiotech, LLC.